



## MUG DISKS

- For in vitro use only -

Catalogue No. DM95

Our MUG Disks are used to aid in the detection of *Escherichia coli*, and to help differentiate between *Streptococcus equisimilis* and *Streptococcus anginosus* (also known as *S. milleri*).

Organisms possessing the enzyme  $\beta$ -glucuronidase can hydrolyze MUG (4-methylumbelliferyl- $\beta$ -D-glucuronide) releasing the fluorescent end product, 4-methylumbelliferone.

Detection of this end product can be made visually under a long-wave (365 nm) UV light source, as 4-methylumbelliferone fluoresces bluish-white to blue-green. Studies conducted by Feng and Hartman showed  $\beta$ -glucuronidase activity in 96% of the strains of *E. coli*, 100% of enterotoxigenic strains of *E. coli*, 17% of *Salmonella* species, and 40% of *Shigella* species.

### Recommended Procedure

1. Place disk in an empty petri dish.
2. Moisten slightly with sterile saline.
3. Using a sterile inoculating loop, inoculate the moistened disk with the organism to be tested.
4. Incubate for 30 minutes at 35°C.
5. Examine under a long-wave UV light for fluorescence.

### Interpretation of Results

Positive: Fluorescence

Negative: No fluorescence

### Quality Control

<u>Organism</u>	<u>Expected Results</u>	
<i>Escherichia coli</i> ATCC 25922	+ve	Fluorescence
<i>Proteus mirabilis</i> ATCC 12453	-ve	No fluorescence

### Storage and Shelf Life

MUG Disks should be stored at -20°C. At this temperature they have a shelf life of 32 weeks from the date of manufacture.

### References

1. Lawrence J, Yajko DM, Hudley WK. J Clin Micro 1985; 22:772-777.
2. Feng PCS, Hartman PA. Appl Environ Micro 1982; 43:1320-1329.

Original: August 2000

Revised / Revisited: October 2014