



PYR DISKS

- For in vitro use only -

Catalogue No. DP95 / DP95K

Our PYR Disks are intended for the preliminary rapid characterization of β -hemolytic, gram-positive cocci.

Organisms possessing the enzyme L-pyrrolidonyl arylamidase can hydrolyze the disk substrate L-pyrrolidonyl- β -naphthylamide to release L-pyrrolidone and β -naphthylamide. Visual detection can be achieved by the addition of PYR / LAP Reagent (Dalynn RP95); the active ingredient, p-dimethylaminocinnamaldehyde, combines with the end-product β -naphthylamide to form a **red** Schiff base. A yellow to red color change indicates a positive reaction.

PYR disks are mainly used for the presumptive identification of group A β -hemolytic streptococci (i.e. *Streptococcus pyogenes*). Yajko et al. reported a high sensitivity (96.4%) and higher specificity (98.7%) for the PYR hydrolysis test as compared to traditional bacitracin susceptibility tests. Among the commonly encountered β -hemolytic streptococci only *S. pyogenes* produces a positive reaction. The accuracy and rapidity of the PYR test allows for quicker reporting of results for those cultures screened for group A streptococci and can result in better patient care and treatment.

PYR disks can also be used for the presumptive identification of *Enterococcus* species. The most accurate presumptive identification of a catalase-negative gram-positive coccus as an *Enterococcus* strain is by demonstrating that the unknown strain is PYR and LAP positive and grows in 6.5% NaCl, and at 45°C.

Recommended Procedure

1. Place a PYR Disk into a suitable sterile container and allow the disk to warm to room temperature.
2. Rehydrate the disk with a small drop of sterile purified water or phosphate buffer (pH 7.5).

3. Inoculate the disk surface with several pure colonies of the unknown test organism derived from an 18-24 hour culture on a blood agar plate.
4. Incubate at room temperature for 10 minutes.
5. Add one drop of PYR / LAP Reagent to the inoculated disk.
6. Interpret results within one minute.

Interpretation of Results

Positive: Cherry red color
Group A β -streptococci (*S. pyogenes*)
and group D enterococci

Negative: Orange or yellow (no change)
Group B streptococci, viridans
streptococci, β -streptococci
(not A, B, or D)

Organism	Catalase	LAP	PYR	Esculin	6.5% NaCl	Vanco
<i>β-Streptococcus</i>						
<i>S. pneumoniae</i>	-	+	V	V	-	S
<i>S. pyogenes</i> (A)	-	+	+	+	-	S
Other β -Strep	-	V	-	-	-	S
<i>Enterococcus</i> spp.	-	+	+	+	V+	S
<i>Abitrophia</i>						
<i>A. adiacens</i>	-	+	V	-	-	S
<i>A. defectiva</i>	-	+	V	-	-	S
<i>Aerococcus</i>						
<i>A. viridans</i>	-	-	+	V	+	S
<i>A. urinae</i>	-	+	-	V	+	S
<i>Alliococcus otitis</i>	+	+	+	V	+a	S
<i>Gemella</i>						
<i>G. hemolysans</i>	-	V	+	-	-	S
<i>G. morbillorum</i>	-	+	+b	-	-	S
<i>Helcococcus kunzii</i>	-	-	+	+	V	S
<i>Lactococcus</i> spp.	-	+	+	+	-	S
<i>Leuconostoc</i> spp.	-	-	-	-	V	R
<i>Pediococcus</i> spp.	-	+	-	+	V	R
<i>Tetragenacoccus</i> spp.	-	+	-	NR	NR	S
<i>Vagacoccus</i> spp.	-	+	+	+	+	S

V = Variable; NR = No Result; S = Susceptible; R = resistant

a = May require 2 to 7 days; b = Weakly positive with large inoculum

- Ensure that the test organism is β -hemolytic, catalase negative, and a gram-positive coccus before performing the PYR test
- *Streptococcus porcinus* and *S. iniae* are animal associated, β -hemolytic species that possess the necessary enzyme to produce a positive reaction
- Enterococci and group A streptococci are both PYR-positive, but differences in colony size and morphology should allow for differentiation
- Ensure that an adequate inoculum is used or false negatives may occur

Quality Control

Organism	Expected Results	
<i>Streptococcus pyogenes</i> ATCC 19615	+ve	Red color change
<i>Streptococcus agalactiae</i> ATCC 27956	-ve	No color change

Storage and Shelf Life

PYR Disks should be stored at 2°C to 8°C, and protected from light. Under these conditions the disks have a shelf life of 52 weeks from the date of manufacture.

Ordering Information

Cat# DP95-25 PYR Disks (25/vial)
 DP95-50 PYR Disks (50/vial)
 DP95K PYR Kit (25 Disks & 3-mL Reagent)

References

1. Godsey J, Schulman R, Enriquez I. The hydrolysis of L-pyrrolidonyl- β -naphthylamide as an aid in the rapid identification of *Streptococcus pyogenes*, *S. avium*, and group D enterococci. Abstr Annu Meet ASM 1981; C84:276.
2. Facklam RR, Thacker LG, Fox B, Enriquez L. Presumptive identification of streptococci with a new test system. J Clin Micro 1982; 15:987-90.
3. Wasilauskas BL, Hampton KD. Evaluation of the Strep-A-Fluor identification method for strep A streptococci. J Clin Micro 1984; 20:1205-1206.
4. Yajko DM, Lawrence J, Nassos P, Young J, Hadley KW. Clinical trial comparing bacitracin with Strep-A-Chek for accuracy and turnaround time in the presumptive identification of *Streptococcus pyogenes*. J Clin Micro 1986; 24:431-4.
5. Gordon LP, Damm MAS, Anderson JD. Rapid presumptive identification of streptococci directly from blood cultures by serologic tests and the L-pyrrolidonyl- β -naphthylamide reaction. J Clin Micro 1987; 25:238-41.
6. Murray PR, Baron EJ, Pfaller MA, Tenover FC, Tenover RH, Eds. Manual of clinical microbiology. 7th ed. Washington, DC: ASM Press, 1999.
7. MacFaddin, JF. Biochemical Tests for the Identification of Medical Bacteria, 3rd ed. Philadelphia: Lippincott Williams & Wilkins, 2000.

Original: January 2002

Revised / Reviewed: October 2014