# **NAD Supplement**

- For in vitro use only -

## Catalogue No.

VN15-05 NAD Supplement [5-mL / 500-mL] (Liquid) LN15-05 NAD Supplement [5-mL / 500-mL] (Lyophilized)

A nutritional supplement used in the production of Haemophilus Test Medium, which is used in the antimicrobial disk diffusion susceptibility testing of *Haemophilus* species. The active ingredient,  $\beta$ -nicotinamide adenine dinucleotide (NAD), is a key nutritional element required by most *Haemophilus* species for growth.

## Active Ingredients per 5-mL Vial (Each vial prepares 500-mL of media)

NAD......7.5 mg

### **Appropriate Commercial Bases**

Manufacturer	Description	Catalogue No.
Difco	Mueller Hinton Agar	225250
BBL	Mueller Hinton II Agar	211438
Oxoid	Mueller Hinton Agar	CM337

### **Reconstitution Procedure**

The lyophilized supplement must be reconstituted prior to use:

- 1. Aseptically add 5.0-mL of sterile purified water to the vial.
- 2. Swirl vial gently until supplement is completely dissolved.

## **Recommended Method for Media Preparation** (**Rose Bengal Agar**)

- 1. Allow one 5.0-mL vial of NAD Supplement to adjust to room temperature prior to using.
- 2. Prepare 500-mL of Mueller Hinton Agar as outlined by manufacturers instructions along with yeast extract (5.0-g/L) and solubilized hematin (0.015-g/L). The medium should be mixed with heat until boiled and transferred to an autoclave for sterilization.
- 3. Autoclave medium for 15 minutes at 121°C.
- 4. After sterilization, cool medium to 45 to 50°C in a warm water bath.
- 5. Mix supplement well before adding. Add the 5.0-mL vial of NAD Supplement to the medium; mix well to incorporate the supplement thoroughly into the medium.
- 6. Pour the completed medium into sterile petri dishes.
- 7. Allow medium to set on a level surface.



## **Quality Control**

The following organisms are used to determine the performance of the completed medium. Inoculate and incubate in a CO<sub>2</sub>-enriched atmosphere at 35°C for 18 to 24 hours.

Organism	Expected Results
Haemophilus influenzae ATCC 49247	Growth
Haemophilus influenzae ATCC 49766	Growth
Haemophilus influenzae ATCC 10211	Growth

#### **Storage and Shelf Life**

Our frozen NAD Supplement has a shelf life of 32 weeks from the date of manufacture when stored at -20°C. The lyophilized supplement has a shelf life of 104 weeks (2 years) from the date of manufacture when stored at 4 to  $8^{\circ}$ C.

Original: December 2006 Updated: December 2006 Checked: January 2014