



OFPBL AGAR

- For in vitro use only -

Catalogue No. PO20

Our OFPBL Agar is used for the selective isolation of *Burkholderia cepacia* from clinical specimens as well as non-clinical samples.

Burkholderia cepacia is an opportunistic pathogen that has been associated with nosocomial infections caused by contaminated equipment, medications, and disinfectants, although the most at risk group are CF patients. Patients with cystic fibrosis have a predisposition for infection and infected patients, if untreated, show a rapid decline in lung function, frequent bacteremia, and death due to lung failure.

OFPBL stands for Oxidation Fermentation basal medium with Polymixin Bacitracin and Lactose. This medium is a modified version of OF basal medium developed by a group of researchers (Welch et al.) in 1987. They demonstrated that OFPBL agar resulted in improved recovery and isolation of *Burkholderia cepacia* when compared to other mediums such as MacConkey agar, XLD agar and various blood agars.

OFPBL contains the nutritional components pancreatic digest of casein and lactose. Lactose is a carbohydrate readily utilized by *Burkholderia cepacia*. The fermentation of lactose results in the release of acid end-products detected by the pH indicator, bromthymol blue, present in the medium. When sufficient acid is produced the medium changes from green to yellow providing the colonies their yellow coloration. The selectivity of the medium owes itself to the presence of the antibiotics polymixin B and bacitracin; together these antibiotics provide good suppression of the bacterial flora present in respiratory secretions and sputum.

Formulation per Litre of Medium

Pancreatic Digest of Casein	2.0 g
Sodium Chloride	5.0 g
Lactose	10.0 g
Dipotassium phosphate	0.3 g
Bromthymol Blue	0.08 g
Agar	15.0 g
Polymixin B	600,000 IU
Bacitracin	200 IU

pH 6.8 ± 0.2

Recommended Procedure

1. Allow plates to adjust to room temperature prior to inoculation.
2. Appropriate patient samples for testing include sputum, bronchial washings, and pharyngeal swabs.
3. Using a direct or diluted inoculum from the sample, perform a four-quadrant streak to obtain well-isolated colonies. If the sample is contained on a swab, roll the swab over a small area near the edge of the plate and proceed to streak for isolation using a sterile loop.
4. Incubate plates aerobically at 30 to 35°C.
5. Examine plates daily for up to 5 days.

Interpretation of Results

Typically, *Burkholderia cepacia* colonies appear as yellow colonies with yellow halos. Most strains will grow in 48 hours but some strains may require up to 5 days for the color development. Incubate plates a minimum of 5 days before discarding.

